

Flexsafe® 2D Bags


Sartorius Flexsafe® 2D bags are known well in the industry for its robustness, strength, and flexibility, making them a leader in the storage and shipment of critical liquids in the biopharma world. From Media/Buffer preparation to Downstream processing, Flexsafe® 2D Bags offer solutions throughout the entire process.



Flexsafe® 2D Bags allow sterile connections in:

- Media Preparation
- Buffer Preparation
- Cell Harvest & Downstream Processing





- USP Class VI
- USP <87> & ISO 10993
- USP <88> & ISO 10993
- Visual Inspection of film, bag, & seal performed by manufacturer
- Endotoxin USP <85> & E.P. 2.6.14
- TSE/BSE—EMA/410/01 & E.P. 5.5.8
- Leak-proof testing performed by manufacturer


Customization Available

Cut down labor hours buy customizing your bag to fit your projects needs. Flexsafe bags come with 3 ports, enabling many different combinations for even the most complex jobs. Build your bag with filters, sterile connectors, sterile disconnectors, and much more. See below for some of Bioprocess Supplies Standardized bags to help get you started.

Contact Bioprocess Supplies to begin designing your customized Flexsafe Bag today!

267-313-4534

INFO@BPSSU.COM



Bag Volume	Part Number	P1	P2	P3
3L	BP-3L-MPC-FS-4 BP-3L-MPC-FS BP-3L-AQG-FS	Sartopore 2 -0.45/0.2um - MidiCaps size 4 3/8" MPC Female and MPC Plug 3/8" AQG	TPE Needleless Injection Port and Cap	3/8" MPC Male and MPC Cap
5L	BP-5L-MPC-FS-4 BP-5L-MPC-FS BP-5L-AQG-FS	Sartopore 2 -- 0.45/0.2 um - MidiCap size 4 3/8" MPC Female and MPC Plug 3/8" AQG		
10L	BP-10L-MPC-FS-4 BP-10L-MPC-FS BP-10L-AQG-FS	Sartopore 2 -- 0.45/0.2 um - MidiCap size 4 3/8" MPC Female and MPC Plug 3/8" AQG		
20L	BP-20L-MPC-FS-7 BP-20L-MPC-FS BP-20L-AQG-FS	Sartopore 2 -- 0.45/0.2um - MidiCap size 7 3/8" MPC Female and MPC Plug 3/8" AQG		
50L	BP-50L-MPC-FS-7 BP-50L-MPC-FS BP-50L-AQG-FS	Sartopore 2 -- 0.45/0.2um - MidiCap size 7 3/8" MPC Female and MPC Plug 3/8" AQG		